

D2-Lampe Shimadzu UV-Serie

Artikelnummer 092600



Produktbeschreibung

D2-Lampe für Shimadzu UV-120, UV-150(-2), UV-160(A), UV-240, UV-250, UV-260, UV-365, UV-625, UV-1200 Serie, UV-1600 Serie, UV-1700 Serie, UV-1800 Serie, UV-2100, UV-2401PC, UV-2501PC, UV-2600, UV-2700, UV-3000, UV-3100, UV-3600, UV-3700, AA-670, AA-6200, AA-6300, AA-6500, AA-6800, CS-920, CS-930, CS-9000

Methode: Photometrie

Vergleichbar mit: 062-65055-05, 062-65062-01, 062-65063-01, 200-75528-SE, 200-75503-01

Brenndauer: 1000 h

Lampen-Typ: Deuterium-Lampe

Produkteigenschaften

Geräte-Hersteller Shimadzu

Geräte-Modell Shimadzu UV Series

Weitere Bilder

